ADVANCED WOUND DRESSING SELECTION GUIDE

ERATE DRAINAGE IMAL DRAINAGE VY DRAINAGE

**Note:** Per patient monthly maximums:

120 Collagen Dressings 120 Foam Dressings

120 Hydrogel Dressings



REV.6.20.24

	9	≡	$\overline{\Box}$				
DRESSING TYPE	NONE	MINI	MODI	HEAV	THICKNESS	MAX ALLOWABLE (PER WOUND)	HCPCS
Calcium Alginate			•	•	FULL	<b>30</b> (once daily) Accurate may supply a combined maximum of 30 dressings coded as A6196, A6197	A6196, A6197, A6199*
Collagen		•	•		FULL	30 (once daily)	A6010*, A6021, A6022, A6023
Composite (bordered)			•	•	ANY	12 (three times per week)	A6203, A6204
Compression Bandages		•	•	•	ANY	4 (once weekly)	A6448-A6452
Contact Layer	•	•	•	•	ANY	4 (once weekly)	A6206-A6208
Foam (as a primary)			•	•	FULL	12 (three times per week)	A6209, A6210, A6211, A6212, A6213, A6214
Foam (as a secondary)				•	FULL	12 (three times per week) For secondary foams, primary dressing FOC must be 3x/wk or less	A6209, A6210, A6211, A6212, A6213, A6214
Gauze (bordered)	•	•	•	•	ANY	30 (once daily)	A6219, A6220
Gauze (impregnated)	•	•	•	•	ANY	30 (once daily)	A6222, A6223
Gauze (non-impregnated; non-bordered)	•	•	•	•	ANY	60 (two times per day)	A6402
Gauze (non-impregnated; bordered)	•	•	•	•	ANY	30 (once daily)	A6220
Gauze (plain rolled/ stretch/bulkee)	•	•	•	•	ANY	30 (once daily)  FOC is dictated by the FOC of the underlying dressing	A6446
Gauze (packing strip, PHMB rolled; impregnated)	•	•	•	•	ANY	30 (once daily)	A6266*, A6407

REV.6.20.24 ADVANCED WOUND DRESSING SELECTION GUIDE								
DRESSING TYPE				THICKNESS	MAX ALLOWABLE (PER WOUND)	HCPCS		
Multilayer (2-, 4- layer)		•	•	• FULL	4 (once weekly)	A6442-A6455		
Hydrocolloid (all sizes)		•	•	ANY	12 (three times per week)	A6234, A6235, A6237, A6238		
<b>Hydrocolloid</b> (filler - e.g., Triad)		•	•	ANY	6 oz. per wound (three times per week)	A6240*		
Hydrogel (bordered sheets)	•	•		FULL	12 (three times per week)	A6245, A6246, A6247		
Hydrogel (non-bordered sheets, impregnated gauzes)	•	•		FULL	30 (once daily)	A6231, A6232, A6233, A6242, A6243, A6244		
Hydrogel (filler)	•	•		FULL	3 oz. per wound (once daily)	A6248*		
Super Absorbent Dressing			•	• FULL	30 (once daily)	A6196, A6197		
Super Absorbent Dressing (when used w/ alginate)				• FULL	30 (once daily)  Accurate may supply a combined maximum of 30 dressings coded as A6196, A6197 per wound	A6196, A6197		
Specialty Absorptive (bordered - e.g., bordered ABD pad)			•	• FULL	15 (every other day)	A6254		
Specialty Absorptive (non-bordered) (e.g., ABD pad. non-bordered composite, non-adherent pad)			•	• FULL	30 (once daily)	A6251, A6252, A6253		
Transparent Film		•		PARTIAL	12 (three times per week)	A6257, A6258, A6259		
Zinc Paste Impregnated Bandage		•	•	• ANY	4 (once weekly)	A6456		

<sup>\*</sup>May be used as a filler in the presence of undermining/tunneling and/or significant depth and be combined with appropriate primary and secondary dressings (requires proper documentation).